

Appendix F

**Hair 5-Panel plus Expanded Opiates**

Drug Name	Initial Screening Cut-Off Limit	Confirmation Cut-Off Limit
<b>AMPHETAMINES</b> <ul style="list-style-type: none"> <li>• AMPHETAMINE</li> <li>• METHAMPHETAMINE</li> <li>• Ecstasy (MDMA)</li> <li>• MDA</li> </ul>	500 pg/mg hair	500 pg/mg hair 500 pg/mg hair 500 pg/mg hair 50 pg/mg hair
<b>COCAINE</b> <ul style="list-style-type: none"> <li>• Cocaine</li> <li>• Benzoylecgonine</li> <li>• Cocaethylene</li> <li>• Norcocaine</li> </ul>	200 pg/mg hair	200 pg/mg hair 200 pg/mg hair 200 pg/mg hair 20 pg/mg hair
<b>OPIATES</b> <ul style="list-style-type: none"> <li>• CODEINE</li> <li>• MORPHINE</li> <li>• HYDROCODONE</li> <li>• HYDROMORPHONE</li> <li>• Oxycodone</li> <li>• 6MAM (Heroin metabolite)</li> </ul>	200 pg/mg hair	200 pg/mg hair 200 pg/mg hair 200 pg/mg hair 200 pg/mg hair 200 pg/mg hair 200 pg/mg hair
<b>PHENCYCLIDINE</b>	300 pg/mg hair	300 pg/mg hair
<b>MARIJUANA METABOLITE</b>	1 pg/mg hair	0.1 pg/mg hair

**Oral Fluid 10-Panel Test**

Drug Name	Initial Screening Cut-Off Limit	Confirmation Cut-Off Limit
<b>AMPHETAMINES</b> <ul style="list-style-type: none"> <li>• AMPHETAMINE</li> </ul>	100 NG/ML	50 NG/ML
<b>METHAMPHETAMINES</b> <ul style="list-style-type: none"> <li>• METHEAMPHETAMINE</li> <li>• MDEA</li> <li>• MDA</li> <li>• MDM</li> </ul>	40 NG/ML	50 NG/ML 50 NG/ML 50 NG/ML 50 NG/ML
<b>BENZODIAZEPINES</b>	1 NG/ML	5 NG/ML
<b>BARBITURATES</b>	20 NG/ML	60 NG/ML
<b>COCAINE METABOLITE</b>	5 NG/ML	8 NG/ML
<b>OPIATES</b> <ul style="list-style-type: none"> <li>• CODEINE</li> <li>• MORPHINE</li> <li>• HYDROCODONE</li> <li>• 6-MAM (Heroin)</li> </ul>	10 NG/ML	40 NG/ML 40 NG/ML 40 NG/ML 4 NG/ML
<b>OXYCODONE</b>	10 NG/ML	40 NG/ML
<b>PHENCYCLIDINE</b>	1 NG/ML	10 NG/ML
<b>MARIJUANA METABOLITE</b>	1 NG/ML	2 NG/ML
<b>METHADONE</b>	5 NG/ML	15 NG/ML

**11 Panel Expanded Opaite and Amphetamine Urine Drug Test**

Drug Name	Initial Screening Cut-Off Limit	Confirmation Cut-Off Limit
AMPHETAMINES <ul style="list-style-type: none"> <li>• AMPHETAMINE</li> <li>• METHAMPHETAMINE</li> <li>• MDMA</li> <li>• MDA</li> <li>• MDEA</li> </ul>	500 NG/ML	250 NG/ML 250 NG/ML 250 NG/ML 250 NG/ML 250 NG/ML
BENZODIAZEPINES	300 NG/ML	300 NG/ML
BARBITURATES	300 NG/ML	300 NG/ML
COCAINE METABOLITE	150 NG/ML	100 NG/ML
OPIATES <ul style="list-style-type: none"> <li>• CODEINE</li> <li>• MORPHINE</li> <li>• HYDROCODONE</li> <li>• HYDROMORPHONE</li> </ul>	300 NG/ML	300 NG/ML 300 NG/ML 300 NG/ML 300 NG/ML
6-MAM (Monoacetylmorphine)	10 NG/ML	10 NG/ML
OXYCODONE	100 NG/ML	100 NG/ML
PHENCYCLIDINE	25 NG/ML	25 NG/ML
MARIJUANA METABOLITE	50 NG/ML	15 NG/ML
METHADONE	300 NG/ML	300 NG/ML
PROPOXYPHENE	300 NG/ML	300 NG/ML
Alcohol (Breath)	0.02%	0.02%

**Blood 10-Panel Drug Test**

Drug Name	Reference Range	Confirmation Threshold
AMPHETAMINES <ul style="list-style-type: none"> <li>• AMPHETAMINE</li> <li>• METHAMPHETAMINE</li> <li>• MDMA</li> <li>• MDA</li> <li>• MDEA</li> </ul>	50 NG/ML	10 NG/ML
BENZODIAZEPINES	20 NG/ML	10 NG/ML 2 NG/ML for: alprazolam, triazolam, midazolam and clonazepam
BARBITURATES	0.1 UG/ML	0.2 UG/ML
COCAINE METABOLITE	25 NG/ML	10 NG/ML
OPIATES <ul style="list-style-type: none"> <li>• CODEINE</li> <li>• MORPHINE</li> <li>• HYDROCODONE</li> <li>• HYDROMORPHONE</li> </ul>	5 NG/ML	1 NG/ML
6-MAM (Heroin)	5 NG/ML	1 NG/ML
OXYCODONE	5 NG/ML	1 NG/ML
PHENCYCLIDINE	8 NG/ML	5 NG/ML
MARIJUANA METABOLITE	5 NG/ML	1 NG/ML
METHADONE	25 NG/ML	1 NG/ML
PROPOXYPHENE	50 NG/ML	50 NG/ML